

Vivago CARE watch – Fact Sheet

Vivago CARE 8001 watch

Vivago CARE 8005 wandering control watch



Order codes	CARE 8001 watch	CARE 8005 watch
Colour	Code	Code
Vivago CARE watch	EU665880010	EU665880050
Vivago CARE watch with red alarm button	EU621980010	EU621980050

Features	CARE 8001 watch	CARE 8005 watch
Manual alarm	x	x
Light indicator as a sign of an alarm	x	x
Automatic alarms	x	x
Connection lost alarm (15 – 240 min)	x	x
Watch off wrist alarm	x	x
Battery alarms (empty / charging / full)	x	x
Wellbeing measurement (activity, sleep, circadian rhythm)	x	x
Automatic self-diagnostics	x	x
Support for wandering control	-	x

Technical specifications

Dimensions & Weight

Weight	~38 g
Diameter	~42 mm
Thickness	~16,5 mm
Total strap length	~23 cm

Materials

Case / Screws	ABS / Stainless steel
Strap / Buckle	Silicone / Aluminium

Battery

Type	Rechargeable, LiPo
Operation time up to (depending on use)	CARE 8001: 2–4 months CARE 8005: 1–1,5 months
Charging time	About 2 hours
Charger	Charger is not included in the delivery, and needs to be ordered separately.

Radio transmission

Data transfer	CARE 8001: Every 20. seconds CARE 8005: Every 20. seconds, Wandering control every 3.s
Range	Up to 30 meters. Note: Concrete, metal, or other structures may hinder transmission of the signal.
Centre frequency	868.95 MHz
Effective radiated power	Under 10 mW

Operating conditions

Operating/ Storage temperature	+ 10 to + 40 °C / + 0 to + 50 °C
Water proof	IPx7 (not for swimming or hot environment like sauna)
Other	Alarm button needs to be pressed for 2 seconds

Vivago

Vivago Oy

Lars Sonckin kaari 14, 02600 Espoo, Finland
Tel. +358 10 2190 610 | email: info@vivago.com | www.vivago.com